


# Classic Mattress

Delta Medical's Classic Mattress consists of a 3 layered foam offering support from the base layer to superior comfort from the low recovery top layer foam.



## Foam Specifications


- Delta Mattress foam, used in the healthcare industry, has an antimicrobial agent built in which guards against a variety of microbes.
- 3 layered foam, top layer being a low recovery foam.



## Options available:



Custom printing



Custom sizing



Tracking



Product Code:  
CL19386ZE1653

- Mental Health Mattress option also available: • **Product Code: CL19386MH1653**
- Fully welded seams
  - No zips
  - No handles
  - With air vents



All seams are welded ensuring it is impervious to fluid penetration.

## Specifications

Product Code	CL19386Z1653
Length	1930mm
Width	860mm
Depth	150mm
Maximum Weight	160kg
Mattress Weight	12kg
Cover Warranty	3 years
Foam Warranty	7 years

## Cover Specifications

### Top cover

- Vapour Permeable
- Multi directional elasticity
- Waterproof
- Antibacterial
- BS 6807 ignition source 5
- BS 7175 ignition source 5
- Fully welded with zip (4 sides) and waterfall flaps
- 3 year guarantee on fabric if maintained correctly
- Care instructions on base cover

### Base Cover

- Grey vinyl non slip for durability and longevity
- Carry handles

Standard Colour



Carrflex Extreme Deep Sea Green



**DELTA**  
MEDICAL