


# Classic Mattress

Delta Medical's Classic Mattress consists of a 3 layered foam offering support from the base layer to superior comfort from the low recovery top layer foam.



## Foam Specifications


- Delta Mattress foam, used in the healthcare industry, has an antimicrobial agent built in which guards against a variety of microbes.
- 3 layered foam, top layer being a low recovery foam.  
50mm Top layer 52/60  
70mm 2nd layer 39/140  
30mm Base layer 36/200



## Options available:



Custom printing



Custom sizing



Tracking



Product Code:  
CL19386ZE1653

- Mental Health Mattress option also available: • **Product Code: CL19386MH1653**
- Fully welded seams
  - No zips
  - No handles
  - With air vents



All seams are welded ensuring it is impervious to fluid penetration.

## Specifications

|                 |              |
|-----------------|--------------|
| Product Code    | CL19386Z1653 |
| Length          | 1930mm       |
| Width           | 860mm        |
| Depth           | 150mm        |
| Maximum Weight  | 160kg        |
| Mattress Weight | 12kg         |
| Cover Warranty  | 3 years      |
| Foam Warranty   | 7 years      |

## Cover Specifications

### Top cover

- Vapour Permeable
- Multi directional elasticity
- Waterproof
- Antibacterial
- BS 6807 ignition source 5
- BS 7175 ignition source 5
- Fully welded with zip (4 sides) and waterfall flaps
- 3 year guarantee on fabric if maintained correctly
- Care instructions on base cover

### Base Cover

- Grey vinyl non slip for durability and longevity
- Carry handles

Standard Colour



Carrflex Extreme Deep Sea Green



**DELTA**  
MEDICAL