

Delta Face Shields




KEY SPECIFICATION

Product code 10500C – Full face shield


Quantity 1pcs / polybag
50pcs / carton

Sizing Universal
33cm W x 24cm L



Elastic strap for universal fit, ease of use and fast removal.

- Elastic strap for universal fit
- Medical grade foam headband, vented for maximum comfort
- Crystal clear anti-fog protective shield
- TGA approved - ARTG# 114415
- All materials used are latex free
- Full length visor for maximum face and neck protection



Vented foam moulds to your head shape giving maximum comfort for all day wear.

In stock and available now



Delta Medical is a Queensland owned and operated business which has been supplying hospitals and other facilities all over Australia for over 23 years. We are fully accredited with ISO9001:2015 and our products are approved by the TGA where required. Our large warehouse in Loganholme allows us to maintain large stocks of our range of products. Delta Medical is committed to providing its customers with high quality products at an affordable price.

📍 90 Chetwynd Street
LOGANHOLME Q 4129

☎ 1800 033 191

✉ info@paradisegroup.com.au

🌐 www.deltamedical.com.au