

Examination Gown

A disposable non-woven, lint free gown which is breathable, strong and fluid repellent offering protection and comfort. Appropriate for extended wear situations involving light to moderate fluid contact.

KEY SPECIFICATION

Product Code	7036WCBS – Small 7036WCBM – Medium
Sizing	S, M
Quantity	50pcs / carton
Material	SMS 30GSM
ARTG#	153273

STRUCTURE

This gown consists of three layers of non-woven fabric. Outside layers are spunbond polypropylene while the middle layer is meltblown.

SPUNBOND

Spunbond is a single layer non-woven fabric which is broadly defined as a web structure consisting of long continuous strand of filaments offering durability and strength with a cloth like feel.

MELTBLOWN

Meltblown is a web of ultra fine fibres sandwiched between two layers of Spunbond offering high filtration and a barrier to liquids.

FEATURES

- Light weight
- Lint free
- Long sleeves
- Knitted cuffs
- Latex free

AVAILABLE COLOURS: ●

