

# Theatre Briefs



## MATERIAL

### SPUNBOND

Spunbond is a single layer non-woven fabric which is broadly defined as a web structure consisting of long continuous strand of filaments offering durability and strength with a cloth like feel.

### AVAILABLE COLOURS:



## KEY SPECIFICATION

<b>Product code</b>	7017 - White Large 7017XL - White X Large 7017XXL - White XX Large 7017XXXL - White XXX Large 7017XXXXL - White 4X Large
---------------------	--

	7017G - Green Large 7017XLG - Green X Large 7017XXLG - Green XX Large
--	---

<b>Material</b>	SBPP - 34GSM
-----------------	--------------

<b>Quantity</b>	50 pcs / polybag 200pcs / carton
-----------------	-------------------------------------

<b>Sizing</b>	L, XL, XXL, XXXL, XXXXL
---------------	-------------------------

